


DermaPal

# Let's Be Clear

Daavlin Phototherapy Solutions for Psoriasis and Vitiligo



for scalp and spot treatments. DermaPal is safe, affordable and powerful.



# DermaPal

# Daavlin DermaPal • Safe • Affordable • Powerful

# Daavlin delivers a clear solution

for psoriasis, eczema and vitiligo... here's how:

### DermaPal clears your skin

- Treatments last for a few seconds to a few minutes, about three times a week.
- A prescription from your doctor is required for this device.
- Psoriasis patients may see a reduction in symptoms such as redness and flaking in as little as 6-8 treatments. Vitiligo patients may see repigmentation begin within 8 weeks.

### DermaPal is an affordable choice

- Daavlin products are reasonably priced.
- Phototherapy products are frequently covered by insurance including Medicare.
- Insurance assistance is available, call Daavlin and let us do the work for you.

### DermaPal Phototherapy is safe and easy to use

- Safe for most patients, including pregnant women and children, as well as elderly and immuno-compromised patients.
- Plugs into a standard household electrical outlet.
- DermaPal has a self-contained digital timer and doesn't require an inconvenient separate timer.

## Patient Friendly

The astonishingly light-weight treatment wand enables any patient to use the DermaPal without tiring. The durable slide-on comb attachment also doubles as a distance guide. UV protective eyewear and a convenient carrying case are included.

### **Powerful**

Designed with high-intensity lamps, DermaPal clears as effectively as powerful top of the line Daavlin products found in your doctor's office.

# <image>



### Daavlin DermaPal Data

Treatment Area	1′′′ x 4.5″	(2.54 cm x 11.4 cm)
Wand Size	8.25" x 1.75"	(21 cm x 5 cm)
Wand Weight	4 oz	(0.1 kg)
Fixture Size	9″ x 5″	(23 cm x 13 cm)
Total Height	2.75″	(7 cm)
Total Weight:	1 lb 12 oz	(0.8 kg)
Available in 110-120 V, 60 Hz and 220-240 V, 50 or 60 Hz		

All Daavlin phototherapy units are 510k cleared by the FDA.